

Bellavita Auto Bath Lifter

HCPCS: E0240+A9270



FIGURE A




FIGURE B



FIGURE C

- The lightest bathlift on the market at only 20.5 lbs
- The Bellavita lowers to 2.3" (Figure A)
- Seat height reaches 18.8" for deeper tubs
- The backrest reclines to 50° at its lowest position
- Completely tool-free set up makes it easy-to-install and operate When not in use, back easily folds down for space-saving storage (Figure B)
- Watertight hand control floats
- The easy-to-use quick charge lithium-ion battery is accommodated in the hand control
- Two-piece construction makes the lift easy to assemble/disassemble for transport
- Padded seat and backrest provide added comfort
- Washable hygienic covers are sold separately
- Hygienic cutaway and slimline profile for easy handling
- 4 quick-release suction cups for anti-slip bath base
- Headrest available as an option
- Simplify transfers, positioning and rotating when using your bath lift with a Vitaturn Swivel and Transfer Aid
- Suction cups adhere to smooth surfaces only, and not textured surfaces

PRODUCTS

Description

477200252	Bellavita Bath Lift w/ White Cover, 1/EA
-----------	--

SPECIFICATIONS

Description

Backrest Dimensions	25.1" (H) x 12.5" - 14.1" (W)
Backrest Reclines	50 degrees
Baseplate Dimensions	12.2" (W) x 22" (L)
Load Capacity Max.	300 lbs.
Seat Dimensions	13.7" (W) x 16.5" (D) x 22.4" (L)
Shipping Weight	25.3 lbs.
Sitting Height	2.3" - 18.8"
Total Weight (excl. hand control)	20.5 lbs.
Warranty	Limited 3 Year

ACCESSORIES

Description

460900205	Bellavita Side Flap, Replacement, 1/ea, EA
460900252	Bellavita Cover Set, Comfort White, 1/bx, BX
460900312	Bellavita Cover Set, Comfort Blue, 1/ea, BX
412100200	Turning Aid/Transfer aid Vitaturn L, 1/ea, EA

PARTS	Description
460900700	Back Assembly Complete,Bellavita,1 ea, EA
460900800	Seat Assembly Complete,Bellavita,1 ea, EA
460900223	BasePlate,New,Bellavita,1/ea, EA
460900211	Stop Lever Bellavita, 1/ea, EA
460900403	Bellavita Battery Charger, 1/bx, BX
460900213	Bellavita Suction Cup 75mm, 1/ea, EA
460900220	Bellavita Suction Cup Set,Screw in,4/bx, BX
460900600	Bellavita Hand Control w/Battery, 1/ea, EA
460900225	BellavitaClip in/slideinSuction cup,4/pk, PKG